


# STANDARD DRUG CHECKING REAGENTS


For the presumptive identification of the presence of a suspected drug




**REBEL**  
Analytics

	The original reagent. First step when testing MDMA and meth. Also good for identifying 2C-B, Mescaline, and other drugs.		Reacts only with secondary amines. Second step when testing MDMA and meth. Ignore any other colors except blue.		Third step when testing MDMA. Can also identify 5-MAPB, 6-APB and other benzofurans along with some other less common drugs.		Reacts only with primary amines. Good to make sure MDMA does not contain amphetamine, MDA, or 2C-B.		Supplemental. Use alongside others to identify some less common drugs. More reactions available online.		Supplemental. Use alongside others to identify some less common drugs. More reactions available online.		Supplemental. Use alongside others to identify some less common drugs. More reactions available online.	
	<b>MARQUIS</b>	0:00 → 0:30 seconds	<b>A SIMON'S B</b>	INSTANT	<b>FROEHDE</b>	0:00 → 0:30 seconds	<b>A ROBADOPE B</b>	0:00 → 0:30 seconds	<b>LIEBERMANN</b>	0:00 → 0:30 seconds	<b>MANDELIN</b>	0:00 → 0:30 seconds	<b>MECKE</b>	0:00 → 0:30 seconds
MDMA								NO REACTION						
MDA			NO REACTION											
5-APB / 6-APB			NO REACTION											
5-MAPB / 6-MAPB							NO REACTION							
Methamphetamine					NO REACTION		NO REACTION						NO REACTION	
Amphetamine			NO REACTION		NO REACTION								NO REACTION	
2C-B			NO REACTION											
2C-I			NO REACTION		NO REACTION				NO REACTION					
Mescaline			NO REACTION				NO REACTION							
Methylone														
N-Ethyl-Pentylone							NO REACTION		NO REACTION					
Mephedrone	NO REACTION				NO REACTION									
PMA	NO REACTION		NO REACTION										NO REACTION	
PMMA	NO REACTION						NO REACTION						NO REACTION	
Oxycodone			NO REACTION				NO REACTION		NO REACTION					
Heroin			NO REACTION				NO REACTION							
Aspirin			NO REACTION				NO REACTION						NO REACTION	
Sugar			NO REACTION				NO REACTION							

	Use Morris to test cocaine and ketamine. Stir the sample with a toothpick for 30 seconds before observing the final color.	
	<div>A MORRIS B</div>	
	Stir for 30 seconds	
Most Drugs	<div></div>	
Cocaine	<div></div>	
Ketamine	<div></div>	
DCK/2-FDCK	<div></div>	

	Use Ehrlich's and Hofmann to test LSD and other tryptamines. The color change can take up to five minutes.	
	<div>EHRlich'S</div>	<div>HOFMANN</div>
	0:00 → 5:00 minutes	0:00 → 5:00 minutes
LSD	<div></div>	<div></div>
DMT	<div></div>	<div></div>
5-Meo-DMT	<div></div>	<div></div>
Most Drugs	<div>NO REACTION</div>	<div>NO REACTION</div>

	Use Folin to identify the presence of piperazine class drugs.	
	<div>A FOLIN B</div>	
	0:00 → 0:30 seconds	
Most Drugs	<div></div>	
BZP	<div></div>	
TFMPP	<div></div>	